

## Alan-Silk Suture



Easy Retrieval Pack (directions for use)

1 Tear open the foil at the notch in the direction shown



2 Pull the torn portion down as shown till the needle is visible



3 Hold the needle with the needle holder and pull the suture out



### ALAN-SILK

- Sterile surgical suture
- Non absorbable
- Multifilament
- Composed of an organic protein (Fibroin)
- Black braided silk
- Needle made of stainless steel, silicone coated. Eases the passage without damaging the tissue

### Comments:

Class III medical device, single use product


Shelf life: 5 years from production date

Should be stored in a dry and clean place, out of direct sunlight and below a temperature of 30°C

Box size: 62x130x62 mm

### Sizes

| Reference | Needle profile | Needle description | Needle length (mm) | Size | Thread length (cm) |
|-----------|----------------|--------------------|--------------------|------|--------------------|
| SUTSI3    | 3/8            | Reverse cutting    | 19                 | 3/0  | 70                 |
| SUTSI3 ½  | 1/2            | Reverse cutting    | 17.4               | 3/0  | 45                 |
| SUTSI4    | 3/8            | Reverse cutting    | 19                 | 4/0  | 70                 |
| SUTSI4 ½  | 1/2            | Reverse cutting    | 17.4               | 4/0  | 45                 |
| SUTSI5    | 1/2            | Reverse cutting    | 16                 | 5/0  | 70                 |



### Packaging

| Number of pieces / box | Number of boxes / export box |
|------------------------|------------------------------|
| 12                     | 216                          |

We are ISO 13485 certified and our products are CE marked

### ALAN & CO

Avenue André Ernst 17 – 4800 Verviers – Belgium

Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69

www.alan.be – info@alan.be

version 10/2019