

Rhinorolls



COMPOSITION	<p>100% polyamide tubular nose dressing containing an absorbent roll of:</p> <ul style="list-style-type: none"> • 25% cotton • 25% viscose • 50% cellulose
RHINOROLL	<ul style="list-style-type: none"> • High elasticity thanks to the polyamide • One size fits all: +/- 42 cm long, absorbent roll of 8 cm long, Ø 15mm • With ear-loops • Used after nose surgery to absorb secretion emission • Can be used as support for decongestant treatment during common respiratory diseases • Non sterile • Shelf life : 5 years from production date

Packaging

References	Description	Pieces/Box	Weight per box	Boxes/Export box	Export boxes/Pallet
RRA20	Rhinorolls	20	+/- 120g	100	10

Comments :

Class I medical device, single use product
Packed in cardboard boxes 23x9x6 cm to protect from dust
Export box sizes: 59x38x72 cm
General tolerance: +/- 10%
Should be stored in a clean and dry place, out of direct sunlight

We are ISO 13485 certified and our products are CE marked