

Alan Lact Fast suture



Easy Retrieval Pack (directions for use)

1 Tear open the foil at the notch in the direction shown



2 Pull the torn portion down as shown till the needle is visible



3 Hold the needle with the needle holder and pull the suture out



ALAN LACT FAST

- Sterile surgical suture
- Absorbable
- Multifilament
- Composed of : 90 % Glycolide, 10 % Lactide : Polyglactin 910
- Violet coloured with D & C Violet No.2
- Coated with combination of Polyglactin 370 and calcium stearate
- Total loss of tensile strength: 10 to 14 days
- Complete absorption after +/- 42 days
- Needle made of stainless steel, silicone coated. Eases the passage without damaging the tissue

Comments:

Class III medical device, single use product


Shelf life: 5 years from production date

Should be stored in a clean and dry place, out of direct sunlight and below a temperature of 30°C

Box size: 62x130x62 mm

Sizes

Reference	Needle profile	Needle description	Needle length (mm)	Size	Thread length (cm)
SUTLF3	3/8	Reverse cutting	19	3/0	70
SUTLF4	3/8	Reverse cutting	19	4/0	70
SUTLF5	1/2	Reverse cutting	16	5/0	70



Packaging

Number of pieces/box	Number of boxes/export box
12	216

We are ISO 13485 certified and our products are CE marked

ALAN & CO

Avenue André Ernst 17 – 4800 Verviers – Belgium

Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69

www.alan.be – info@alan.be

version 10/2019