

## Astringent solution



<b>COMPOSITION</b>	<ul style="list-style-type: none"> <li>• Octadecahydrated aluminium sulphate, PH.Eur. (25%)</li> <li>• Distilled Water (75%)</li> <li>• <b>Without epinephrine</b> : no risk of circulatory disorders</li> </ul>
<b>ASTRINGENT SOLUTION</b>	<ul style="list-style-type: none"> <li>• Astringent and haemostatic solution used for treatment of gingival retraction together with our retraction cord, knitted or braided</li> <li>• For dental use only</li> <li>• Less sensitive to hydrolysis than Aluminium Chloride (more stable)</li> <li>• Colourless, non-staining</li> <li>• Shelf life : 3 years from production date</li> </ul>

### Packaging :

Reference	Description	Solution/Bottle	Bottles/Case	Cases/Box
ALST20	Glass bottle	20 ml	1	49

### Comments:

Cannot be used in case of severe ulcerous gingivitis

Precaution: slightly acidic salt, rinse thoroughly in case of contact with eyes

Sold in glass bottle in order to keep the product stable

Weight of the bottle: 41,5 g

Case sizes: 30 x 30 x 63 mm

Export box sizes: 225 x 225 x 80 mm

General tolerance: +/- 10 %

Should be stored in a clean and dry place, out of direct sunlight and at room temperature

**We are ISO 13485 certified and our products are CE marked**