

## Podorolls



<b>COMPOSITION</b>	<ul style="list-style-type: none"> <li>• 100% absorbent cotton</li> <li>• H2O2 bleached (non chlorine bleached, free of any optical brightener)</li> </ul>
<b>PODOROLLS</b>	<ul style="list-style-type: none"> <li>• Stability of form (high cohesion)</li> <li>• Do not slip</li> <li>• Hygienic because it's a disposable product</li> <li>• Fit closely the toe</li> <li>• For all chiropodist works, an alternative to the rigid toe separators</li> <li>• Can be used to remove excess product</li> <li>• Shelf life : 5 years from fabrication time</li> </ul>

### Sizes and packaging

Reference	Description	Rolls/bag	Pieces /bag	Bags /boxes	Boxes /Pallet
CRN1/15 CM	Podorolls 50 pcs, length 15 cm	2	100	30	28

### Comments :

Class I medical device  
 Packed in resealable polyethylene bags to protect from dust and humidity.  
 Other packaging available on request  
 Export box size : 59x39x24 cm  
 General tolerance : +/- 10%  
 Should be kept in a dry and clean place

**We are ISO 13485 certified and our products are CE marked**

