TECHNICAL DATA SHEET

dental & medical	Alan Lact Fast suture					
MATERIAL	 Composed of: 90% Glycolide, 10% Polyglactin 910 Violet coloured with D & C Violet No.2 Coated with combination of Polyglactin 370 and calcium stearate Needle made of stainless steel, silicone coated 					
PRODUCT INFORMATION	 Sterile surgical suture Absorbable Multifilament Total loss of tensile strength after 10 to 14 days Complete absorption after +/- 42 days Needle eases the passage without damaging the tissue Shelf life: 5 years from production date Class III medical device sterile, single use product General tolerance: +/- 10 % 					
PACKAGING	 12 sterile pouch/box to protect from dust Alan brand only Box dimensions: 62x130x62 mm Export box dimensions: 60x40x60 cm, 216 boites/carton Should be stored in a clean and dry place, out of direct sunlight and below a temperature of 30°C 					

STANDARDS PACKAGING

Reference	Needle profile	Needle description	Needle length	Size	Thread length	Ouverture facile	
SUTLF3	3/8	Reverse cutting	19 mm	3/0	70 cm	Déchirer la feuille à	TEAR
SUTLF4	3/8	Reverse cutting	19 mm	4/0	70 cm	l'entaille dans le sens	1351
						indiqué	212
	1/2 CIRCLE	3/8 CIRCLE		V		Abaisser la partie déchirée jusqu'à œ que l'aiguille soit visible	7 the
		$ \rightarrow $	REVERSE CUTTING A triangular body offers on the outer convex cur		edge	Utiliser le porte aiguille et retirer la suture	2 ym

ALAN & CO Avenue André Ernst 17 – 4800 Verviers – Belgium Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69

We are ISO 13485 certified and our products are CE marked